



Invacare®

Comet® Alpine+

Challenge the outdoors with power,
control and stability



For those who enjoy adventures across terrains too challenging for a traditional scooter, the **Comet⁺ Alpine⁺** offers the perfect solution. Its robust and powerful driving characteristics gives you the freedom to discover the great outdoors with minimal effort in continuous comfort and safety.

As part of the **Comet** family, the **Comet⁺ Alpine⁺** comes packed with all of its well-known safety features and customisation capabilities, offering a reliable, powerful and stylish scooter.

- Powerful and robust
- Excellent outdoor capabilities
- Safety at its core

Powerful and robust

Whether you're tackling uneven terrain or steep and challenging hills, the **Comet⁺ Alpine⁺** will offer excellent maneuverability and power without compromising safety and style. Thanks to its powerful motors and superior suspension system, the **Comet⁺ Alpine⁺** lets you challenge the great outdoors with minimal effort and complete control.

Features and options



Seat suspension module

Offers additional comfort by adding an optional seat suspension.



Seat riser

Electronically increase the seat height to achieve a preferred driving height.



Off-board charging

An optional off board-charger allows you to conveniently charge the batteries away from the scooter.

Excellent outdoor capabilities

When faced with wet or uneven ground, you can rest assured that vital components will be protected from water and dirt thanks to the in-built anti splash guards. Feel stable and at ease when driving at high speeds with directional stability to keep you moving in a straight line and automatic speed reduction to slow you down when driving round corners.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Comet⁺ Alpine⁺ 510 mm	Comfort and Deluxe: 470 mm Premium: 460 mm	Min. 440 mm Max. 510 mm	Comfort: 500 mm Deluxe: 490 mm Premium: 630 mm	665 - 685 mm	1450 mm	143 kg	approx. 63 kg
Comet⁺ Alpine⁺ 160 kg	2 x 12 V/75 Ah (C20) AGM	2750 mm	100 mm	100 mm	10°	49 km	10 km/h
	Colour of shrouds				Upholstery		
Comet⁺ Alpine⁺ 10km/h: 650 W / 1800 W ¹	Marble White Sapphire Blue Jade Green Luxor Beige Crystal Copper Carmine Red Iconic Chrome Diamond Silver Electric Blue Jasper Red				Black Vinyl		

1. Peak performance

Invacare Australia Pty Ltd:
ABN 45 074 676 378
PO Box 5002, 1 Lenton Place, North
Rocks, NSW 2151, Australia.
Freephone: 1800 460 460,
Fax: 02 8839 5353
Email: orders@invacare.com.au
Web: www.invacare.com.au

Invacare New Zealand:
PO Box 62-124, 4 Westfield
Place, Mt. Wellington,
Auckland, New Zealand.
Freephone 0800 468 222,
Freefax 0800 807 788
Email: sales@invacare.co.nz
Web: www.invacare.co.nz

Invacare Asia Ltd:
PO Box 5002, 1 Lenton Place,
North Rocks, NSW 2151,
Australia.
Phone: +61 2 8839 5330,
Fax: +61 2 8839 5331
Email: asiasales@invacare.com
Web: www.invacare.com

Trusted Solutions,
Passionate People

