

Introducing CairMax by Novis  
– a non-powered hybrid foam air  
combination pressure redistribution surface. It  
is suitable for the prevention and treatment of skin  
breakdown and pressure injuries in patients at risk.

### General

	Standard	King Single
<b>CODE</b>	FAMCM-R02	FAMCM-R02K
<b>CAPACITY (SAFE WORKING LOAD)</b>	230 kg	
<b>WARRANTY</b>	2 year	
<b>PRESSURE INJURY RISK</b>	At risk to high risk	
<b>ARTG</b>	292670	
<b>CARTON DIMENSIONS (L x W x H) in mm</b>	2270 x 900 x 180 mm	2270 x 1070 x 180 mm
<b>DIMENSIONS (L x W x H) in mm</b>	2000 x 850 x 200 mm	2000 x 1050 x 200 mm
<b>MATTRESS WEIGHT</b>	15 kg	20.4 kg
<b>CELLS</b>	<b>NUMBER</b>	10
	<b>FOOT ZONE</b>	Tapered fFoam
	<b>CELL MATERIAL</b>	TPU laminated nylon
	<b>CELL FOAM INSERT MATERIAL</b>	High resilience PU
	<b>MATTRESS AIR SUPPLY LINE</b>	Heavy duty PVC piping
<b>TOP COVER</b>	<b>MATERIAL / TYPE</b>	PU-laminated nylon
	<b>FIRE RETARDANCY</b>	BS7175 Crib 5
	<b>MAXIMUM WASH TEMPERATURE</b>	95° C
<b>BASE</b>	<b>ATTACHMENT</b>	Dual zippers
	<b>SEAM CONSTRUCTION</b>	External seams RF welded
<b>BASE</b>	<b>MATERIAL / TYPE</b>	PVC-laminated polyester

### Support Surface



THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.