

Introducing ProCair
by Novis – alternating
solutions for the prevention and
treatment of pressure injuries and other
preventable complications of immobility.

		Standard size	King Single Size	
General	CODE	APMPC-R01	APMPC-R01K	
	CAPACITY (SAFE WORKING LOAD)	40 - 200 kg		
	SYSTEM ALTERNATION	1 in 2		
	WARRANTY	3 year		
	CARRY BAG	Deluxe carry bag with padded pump and mattress compartments		
	PRESSURE INJURY RISK	High to very high risk		
	COMPLIANCE	IEC 60601-1, IEC60601-1-2, IEC60601-1-11		
	ARTG	289458		
	CARTON DIMENSIONS (L x W x H) in mm	1120 x 430 x 400 mm		
	DIMENSIONS (L x W x H) in mm	2000 x 880 x 200 mm	2000 x 1050 x 200 mm	
MATTRESS WEIGHT	8 kg	8.5 kg		
CPR TYPE	QuickTwist CPR at head end			
Support Surface	CELLS	NUMBER	18	
		STATIC HEAD CELLS / ZONE	3	
	CELLS	CELL-IN-CELL	Yes	
		CELL MATERIAL	TPU-laminated nylon	
		MATTRESS AIR SUPPLY LINE	Welded TPU air manifold and CPC quick release connectors	
	TOP COVER	MATERIAL / TYPE	Dark blue, polyurethane-laminated nylon – water resistant, vapour permeable multi-stretch	
		FIRE RETARDANCY RATE	BS7175 Crib 5	
		MAXIMUM WASH TEMPERATURE	95° C	
ATTACHEMENT		Dual zippers		
	SEAM CONSTRUCTION	RF welded on all external seams		
BASE	MATERIAL / TYPE	Black PVC laminated polyester		
	HINGE	No		
	FIXING STRAPS	Webbing at 4 corners, 2 at middle		

Control
Unit

DIMENSIONS (L x W x H) in mm	350 x 135 x 224 mm	
WEIGHT	3.1kg	
POWER	110 - 220V 50/60Hz	
COMPRESSORS	Single	
OPERATION	Full digital control	
CYCLE TIME	12 minutes	
NOISE LEVEL, WITHIN 2 METERS	32 db(A) ± 1	
MIN / MAX PRESSURE	20-60mm HG +/- 6mm HG	
ELECTRICAL PROTECTION TYPE	Class II Type BF	
OPERATING ENVIRONMENT	Operation 15°C to 35°C Storage 5°C to 60°C	
PRESSURE SETTING	AUTOMATIC	Yes – cyclical
	MANUAL	5 levels
FOWLER FUNCTION	Yes, manual	
BACKUP BATTERY	No	
ALARMS	Low Pressure High Pressure Alternation Failure Startup Failure Power Failure	



THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.