

The Koval Mover is an excellent substitute for a wheelchair with its rotating padded seats and a knee support, which allow comfortable sitting.

The Mover allows fast, stable and safe transport of a patient to the bed, chair or the toilet and back. The rear castors can be fixed by means of a brake.

The Mover is useful for transportation of patients, who have difficulties walking or standing.

Mover Flexible Base Features

- The maximum dimensions of the walking trolley are as follows:
104 cm (height), 97 cm (length),
58.5 cm (width)
- Castors with a diameter of 75 mm,
rear castors are equipped with
a brake
- Trolley weight: 19.6 kg
- Safe working load: 150 kg

Rotating padded seats, excellent substitute for a wheelchair.



Mover EL Features

- Combines features of Mover with
stand-up lifter function
- The maximum dimensions of the
walking trolley are as follows:
104 cm (height), 97 cm (length),
58.5 cm (width)
- Castors with a diameter of 75 mm,
rear castors are equipped with
a brake
- Trolley weight: 19.6 kg
- Safe working load: 150 kg

