

Shower Trolley L1

Features

- Shower trolley for transportation with a one-piece stretcher top
- The maximum dimensions of the shower trolley are as follows: 207 cm (length) and 81 cm (width)
- Stretcher top length: 190 cm, stretcher top width: 65 cm
- Fixed 3° stretcher top inclination
- Foldable side rails (height: 20 cm)
- Hydraulically or electrically lowering/raising of the stretcher top via a console (lowering/raising range: 54 cm to 97 cm)
- Castors with a diameter of 150 mm
- Central braking system
- Easy steering system (for ease of maneuvering the trolley's front left castor can be locked in the direction of steering)
- Metal base cover
- Water-resistant shower mattress (height: 27 cm)
- Mattress integrated drain with a 100 cm long flexible hose
- Water-resistant head pillow – semi-cylinder
- Trolley weight: 86 kg
- Safe working load: 170 kg

One-piece stretcher top.



Options

- Castors with diameter of 200 mm

- Higher water-resistant shower mattress

- 200 cm long flexible hose

- Antistatic castors
