

Ultra Compact Wheeled Walker



The Ultra Compact Wheeled Walker offers all the strength, design and integrity of a regular Wheeled Walker with the added versatility of a compact folding design that allows it to be easily transported or neatly packed away when not in use.

Features

- Compact folding design for easy transportation and storage
- Innovative folding mechanism
- 8" (203mm) all terrain wheels
- Lightweight aluminium frame for strength and stability
- Height adjustable handles with moulded handgrips
- Soft touch braking mechanism with parking brake function
- Large portable shopping bag for added convenience
- Easily removable backrest without the use of tools for easy transport
- Suitable for indoor and outdoor use
- Tilt assist pedals for navigation of obstacles

Specifications

CODE	BL0063
WEIGHT	7.3 kg
OVERALL HEIGHT	785 - 935 mm
SEAT HEIGHT	560 mm
SEAT WIDTH	450 mm
WHEELS	8" (200 mm)
FRAME MATERIAL	Aluminium
SAFE WORKING LOAD	136 Kg
CONFORMS TO	ISO 11199-2:2005
ARTG	176830
WARRANTY	1 Year