

# ALTO Lift™ MkII 200

Heavy Duty Patient Lifting Hoist

ALTO Lift™ MkII  
Now With  
Extended  
Boom



Featuring Linak Jumbo Care advanced diagnostic control system  
Powered by LA44 European Standards compatible linear actuator

**Innovation in  
Patient Handling**

**Allegro**  
CONCEPTS PTY LTD

# ALTO Lift™ MkII 200

**ALTO Lift™ MkII 200:** A medium sized heavy duty mobile patient lifting hoist with a safe working load of 200kg and a compact ergonomic design.

ALTO Lift™ MkII 200, whilst being compact, provides superior patient lift space through the use of SPACE LINK™ technology.

ALTO Lift™ MkII 200 is the lifter of choice for institutional use or home care where a high capacity lift is required.



### SPACE LINK™

- breakthrough for optimized patient space
- enhances floor lift capacity




Standard Yoke




4 Point Yoke



Tri Hook Yoke



Pivot Frame



Powered Pivot Frame



Weigh Scale

### ALTO Lift 200 Specifications

Safe Working Load	200kg
Base Length	1150mm
External Base Width Closed	700mm
External Base Width Open	1130mm
Internal Base Width Closed	600mm
Internal Base Width Open	1030mm
Base Width Adjustment	Electric
Under Bed Clearance	130mm
Minimum Height Yoke	300mm
Maximum Height Yoke	1920mm
Average Lift Range*	1620mm
Lifter Weight	46kg
Minimum Storage Height	1350mm
Under Base Height	55mm
Maximum Hoist Reach	1060mm
Minimum Hoist Reach	610mm
Hoist Reach at Maximum Height	515mm
Mass of Heaviest Part	24kg (Mast assembly)

\* With general purpose large sling.

The Allegro Concepts range of patient lifting equipment including lifters, slings and accessories has been designed and tested to comply with: AS ISO 10535-2002 Hoists for the transfer of disabled persons – Requirements and test methods.

### Allegro Concepts Pty Ltd

30-32 James Street Lidcombe NSW 2141  
 Phone: (02) 9749 7812 Fax (02) 9749 2144  
 Email: sales@allegroconcepts.com.au

**Allegro**  
 CONCEPTS PTY LTD  
[www.allegroconcepts.com.au](http://www.allegroconcepts.com.au)