

High Back Seating

Chairs with features to accommodate the frail user



The BC4 High Back Day Chair

Features:

- Versatile and Stylish Patient Chair
- Full memory foam layer on seat and backrest.
- Luxurious pressure relieving comfort.
- Built-in postural and ergonomic design with fully welded steel frame for strength
- Adjustable Seat Height and Seat Depth
- 10 year structural warranty
- ARTG number: **218697**

Optional

- Optional Timber Arms



Memory Foam Layer to Seat and Back



Optional Timber Arms

Standard Vinyls and Fabrics



Ink Vinyl



Autumn Vinyl



Light Grey Vinyl



Kiwi Vinyl



Latte Vinyl



Hollywood Vinyl



Red Vinyl



Oatmeal Fabric



Charcoal Fabric



Demin Fabric

The BC4 High Back Day Chair Specifications	
Overall Chair Width:	630mm
Seat Width:	510mm
Overall Depth:	780mm
Seat Depth:	430mm to 530mm
Overall Height:	990mm to 1180mm
Seat Height:	440mm to 590mm
Back Height:	720mm
Safe Working Load (SWL):	180kg