

# HEAVY DUTY WALKER

Aluminium

## FEATURES

- Lightweight aluminum frame
- Height adjustable front wheels and rear legs
- Dual button folding system
- Extra wide frame
- Different rear foot options available



## SPECIFICATIONS

Description	Overall Width	Overall Length	Overall Height	Product Weight	Max User Weight
Heavy Duty Walker KA353HD	680mm	560mm	820mm to 1000mm	4.2kg	227kg

Improving your  
life everyday™

[kcare.com.au](http://kcare.com.au)  
1800 783 783

**KCare**  
Healthcare Solutions