

DELUXE SEAT WALKER 6"

Soft touch locking mechanism

Durable construction.
Easily operated.

FEATURES

- Strong aluminum frame
- Light touch hand brake
- Easy fold padded backrest
- Height adjustable handles
- Soft padded seat provides comfort
- With storage bag
- Can be folded in to a compact size for packing
- 6" PVC wheels.



Foldable backrest



Soft Touch Braking System

SPECIFICATIONS

Description	Overall Width	Overall Length	Overall Height	Seat Width	Seat Height	Wheel Diameter	Product Weight	Max User Weight
Deluxe Seat Walker Red KA365RD	620mm	650mm	840mm to 940mm	445mm	550mm	6 inch (152mm)	6.2kg	160kg
Deluxe Seat Walker Blue KA365BD	620mm	650mm	840mm to 940mm	445mm	550mm	6 inch (152mm)	6.2kg	160kg

Improving your
life everyday™

kcare.com.au

KCare
Healthcare Solutions