

# Bowls Walker Deluxe

## FEATURES

- Ultimate safety and strength to support the user
- Extra large seat for ultimate comfort
- Adjustable to suit all users of various size
- Ergonomic design to facilitate support for the user
- Easily foldable frame, lightweight aluminium frame
- Soft Touch Braking system
- Extra soft PU tyres as to not damage grass
- Ideal for bowling, golfing and grass activities
- Privacy pocket under seat



### 48mm Extra Soft Wheels

48mm wide wheels to comply with the Bowls Australia Artificial Device Policy.



### Soft Touch Braking System

Soft touch braking system allows the user to securely lock the wheels by applying very little pressure to the push down brakes. Ideal for constant use to reduce strain when on the green.



### Lightweight & Foldable

Folds down into a convenient size for transportation or storage, and is light enough to be carried when necessary, lifting from one surface to another.

Description	Seat Width	Seat Depth	Wheel Thickness	Wheel Diameter	Seat Height	Min Handle Height	Max Handle Height	Max User Weight
Bowls Walker Deluxe KA367BW	415mm	360mm	48mm	8 inch (203mm)	520mm	780mm	920mm	200kg