

PREMIUM SEAT WALKER

200kg weight capacity

With a Soft Touch locking mechanism, durable construction and easily folded, offering support and safety for the user

FEATURES

- Adjustable to suit all users of various size.
- Ultimate safety and strength to support the user.
- Ergonomic design to facilitate spacious support for the user.
- Lightweight aluminum frame
- Easily foldable frame



SPECIFICATIONS

Description	Overall Width	Overall Length	Overall Height	Seat Width	Seat Depth	Seat Height	Wheel Diameter	Product Weight	Max User Weight
Premium Seat Walker Red Medium KA367MR	620mm	690mm	830mm to 950mm	420mm	350mm	550mm	7 inch (178mm)	8.0kg	200kg

Improving your life everyday™

kcare.com.au

KCare
Healthcare Solutions